

dcanon sintaxis

from /usr/local/bin/dcanon

```
#!/bin/sh
#
# Usage: $0 srcdir dstdir [options [newpatientname newpatientid]]
#
# where options is a SINGLE string, either empty or possibly containing:
#
#     nodesc - remove descriptions (equals
nopatientdesc,nostudydesc,noseriedesc,noimagedesc)
#     nopatientdesc - remove descriptions at patient level
#     nostudydesc - remove descriptions at study/procedure level
#     noseriedesc - remove descriptions at series and procedure step
level
#     noimagedesc - remove descriptions at image level
#     nomanuf - remove manufacturer
#     noretired - remove retired ACR-NEMA elements or other elements that
are undesirable
#     noclinicaltrial - remove clinical trial related attributes
#     keepprivate - do not remove private elements
#     keepsafeprivate - keep safe private attributes only, e.g. Philips
private PET SUV factor
#     keepuids - do not replace UIDs
#     keeppet - keep parameters for PET SUV calculation
(age,sex,size,weight)
#     verify - verify afterwards
#     diff - show difference afterwards
#     mvdatedesc - more results to descriptive hierarchy (default)
#     mvid - more results to hierarchy based on IDs
#     mv8only - more results to 8 character A-Z0-9_ component names
#     nomove - do not move into hierarchy
#     none - dummy option to allow specification of name and id arguments
when no options
#
# e.g. dcanon src dst nodesc,nomanuf Smith^Mary 1234
```

dcm2niix

<https://github.com/rordenlab/dcm2niix>

CODECS chart

https://www.nitrc.org/plugins/mwiki/index.php/dcm2nii:MainPage#DICOM_Transfer_Syntaxes_and_Co

[mpressed_Images](#)

Editar DICOMs

<https://support.dcmtdk.org/docs/dcmtdify.html>

From:

<http://detritus.fundacioace.com/wiki/> - **Detritus Wiki**

Permanent link:

<http://detritus.fundacioace.com/wiki/doku.php?id=neuroimagen:dcanon>

Last update: **2022/01/07 10:09**

